

QuickVue® Dipstick Strep A Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue Dipstick Strep A test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current package insert instructions and/or procedure card included in the kit.

1. How often does the manufacturer recommend that external controls be run on the kit?

- a) Once every 20 tests
- b) Once per each new shipment of kits
- c) Every first sample tested each day
- d) Once per each untrained operator, once for each new shipment of kits-provided that each different lot in a shipment is tested

2. At what temperature do the kits need to be stored?

- a) 2-8°C (refrigerated)
- b) 15-30°C (room temp.)
- c) 2-30°C (refrigerated or room temp.)
- d) -20°C (frozen)

3. For proper test performance, which swab should be used when collecting a throat swab specimen?

- a) The swabs provided in the kit.
- b) Any type of throat swabs
- c) Any Dacron swabs
- d) Cotton tipped swabs

4. What is the read time for this kit?

- a) 5 minutes
- b) 6 minutes
- c) 10 minutes
- d) 20 minutes

5. Can the test be read after the read time?

- a) No
- b) Yes, up to 10 minutes after the read time
- c) Yes, up to 5 minutes after the read time
- d) Yes, any time after the read time

6. How long after adding Reagent A and B together in the tube is it still good?

- a) 8 hours
- b) 48 hours
- c) Immediately after adding both reagents
- d) 24 hours

7. What color should the solution be when Reagents A and B are mixed?

- a) Clear
- b) Any shade of green
- c) Dark green
- d) Yellow

8. What is the incubation time when the swab sample is in the Reagent A and B mixture?

- a) 5 minutes
- b) There isn't an incubation time
- c) 1 minute
- d) 2 minutes

9. The throat swab specimen is good for how long in any clean, dry plastic tube or sleeve before processing?

- a) 24 hours
- b) 8 hours
- c) 48 hours
- d) 72 hours

10. If culture is desired after throat swab is collected, what is the process to follow?

- a) Lightly streak the swab on a 5% sheep blood agar plate before using swab in the QuickVue Dipstick Strep A test or obtain dual throat swab specimens
 - b) There is no process to follow
 - c) Culture is not necessary
-

QuickVue® Dipstick Strep A Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue Dipstick Strep A test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current package insert instructions and/or procedure card included in the kit.

1. How often does the manufacturer recommend that external controls be run on the kit?

- a) Once every 20 tests
- b) Once per each new shipment of kits
- c) Every first sample tested each day
- d) Once per each untrained operator, once for each new shipment of kits-provided that each different lot in a shipment is tested

2. At what temperature do the kits need to be stored?

- a) 2-8°C (refrigerated)
- b) 15-30°C (room temp.)
- c) 2-30°C (refrigerated or room temp.)
- d) -20°C (frozen)

3. For proper test performance, which swab should be used when collecting a throat swab specimen?

- a) The swabs provided in the kit.
- b) Any type of throat swabs
- c) Any Dacron swabs
- d) Cotton tipped swabs

4. What is the read time for this kit?

- a) 5 minutes
- b) 6 minutes
- c) 10 minutes
- d) 20 minutes

5. Can the test be read after the read time?

- a) No
- b) Yes, up to 10 minutes after the read time
- c) Yes, up to 5 minutes after the read time
- d) Yes, any time after the read time

6. How long after adding Reagent A and B together in the tube is it still good?

- a) 8 hours
- b) 48 hours
- c) Immediately after adding both reagents
- d) 24 hours

7. What color should the solution be when Reagents A and B are mixed?

- a) Clear
- b) Any shade of green
- c) Dark green
- d) Yellow

8. What is the incubation time when the swab sample is in the Reagent A and B mixture?

- a) 5 minutes
- b) There isn't an incubation time
- c) 1 minute
- d) 2 minutes

9. The throat swab specimen is good for how long in any clean, dry plastic tube or sleeve before processing?

- a) 24 hours
- b) 8 hours
- c) 48 hours
- d) 72 hours

10. If culture is desired after throat swab is collected, what is the process to follow?

- a) Lightly streak the swab on a 5% sheep blood agar plate before using swab in the QuickVue Dipstick Strep A test or obtain dual throat swab specimens
- b) There is no process to follow
- c) Culture is not necessary



To Order Please Contact:
American Screening Corporation
6658 Youree Dr., Ste 180 PMB 404
"Complete Drug Testing Solutions"
Shreveport, LA 71105
Ph: 318-227-4994 Fax 318-227-4996
www.americanscreeningcorp.com