


OneScreen™

Court Drug Testing



9 Panel in stock
12 panel special order



OneScreen™ Multi-Drug Test Cup

Court Drug Testing Made Simple.

- Lower Drug Level Detection
- Room Temperature Storage
- Detects 3-12 Drugs
- Built-in Validity Test
- 24-Month Shelf Life
- Simple Procedure
- Fast Results
- Patent-Pending Design
- Automatic Split-Specimen With No Urine Contamination



6658 Youree Dr Ste 180 PMB 404 Shreveport, LA 71105


Toll Free: 1-866-526-2873 Ph: 318-798-3306 Fax: 318-798-3386 www.americanscreeningcorp.com

ONE SCREEN™

OneScreen™ Multi-Drug Test Cup

Court drug testing made Simple.

Detect the presence of drugs and adulterants in minutes. It's drug testing efficiency you can count on.



▲
OneScreen 12-Panel Cup
Shown With An Example
Negative Result.

Integrated Testing Solution


Temperature Strip
Two-part Result Form Available

Security Seal
Photocopy Template


Drug Name	Abbreviation	Cutoff Level
• Amphetamine	AMP	300ng/mL
• Barbiturates	BAR	200ng/mL
• Benzodiazepines	BZO	200ng/mL
• Cocaine	COC	100ng/mL
• Marijuana	THC	40ng/mL
• Methadone	MTD	200ng/mL
• Methamphetamine	mAMP	1000ng/mL
• Opiates	OPI	100ng/mL
• Oxycodone	OXY	100ng/mL
• Phencyclidine	PCP	25ng/mL
• Propoxyphene	PPX	300ng/mL

Instructions



1 Collect Specimen




2 Peel Back Label



3 Read the Results



Read results at 5 minutes



TEST	ABNORMAL (LOW)	NORMAL	ABNORMAL (HIGH)
Oxidant/PCC		Negative	
S.G.	1.000	1.003 1.005 1.015 1.025	≥ 1.035
pH	2.0 3.0	4.0 7.0 9.0	10.0 11.0 12.0

SAMPLE ONLY-Color chart NOT Exact