

ALCOMATE™ ACCUCELL (MODEL AL9000)

The AlcoMate AccuCell is our most accurate breathalyzer that provides you with results equivalent to much more expensive law enforcement patrol equipment. The AccuCell combines state-of-the-art fuel-cell alcohol sensing technology with the ease and convenience of one-button operation. Featuring precision and accuracy ratings up to 5 times higher than competing products, the AccuCell is the ideal choice for those seeking the best-of-the-best.

Features

- Fuel-Cell Alcohol Sensor for Maximum Accuracy
- DOT Approved for Law Enforcement (Accuracy +/-0.005)
- Easy One-Button Operation
- Back-Lit 4-Digit (0.000%) LCD Results
- Active Pressure Sensor / Deep-Lung Air Sampling
- Cumulative Test Counter / Audible Alerts / Low Battery Indicator
- One-Way-Flow Mouthpieces
- Self-Diagnostic System / Auto Power Off



Device Dimensions:
4.25" x 2.0" x 0.75"



Retail Case features a sturdy, re-inforced ABS exterior with a thickly padded interior for superior protection of AccuCell and accessories.

MSRP
\$299.95

Retail Case includes: Breathalyzer, User's Manual, 2 x AAA Batteries, 5 x Mouthpieces, Handstrap, Carrying Pouch

Target Users

Government Approvals
for Law Enforcement

DOT (Department of Transportation)
NHTSA (National Highway Traffic &
Safety Administration)



Trusted for Use By

- Law Enforcement
- Medical Professionals
- School Administrators
- Corporate Employers
- Drug Rehabilitators
- Prosumers / General Consumers

6658 Youree Dr Ste 180 PMB 404
Shreveport, La 71105
www.americanscreeningcorp.com
Sales@americanscreeningcorp.com
Toll Free: 1-866-526-2873

Ph: 318-798-3306 Fax: 318-798-3386

Specifications

Alcohol Sensor	Top-Grade Fuel-Cell Sensor
Detection Range	0.000% - 0.400% BAC
Calibration	Required Once Yearly
Breath Sampling Time	4+ seconds (Deep-Lung Air)
Sensor Accuracy	+/- 0.005
Display	4-Digit Decimal LCD (0.000%)
Battery	Two 1.5V AAA alkaline (included)
Weight	1.8 oz (without batteries)
Dimensions	4.25 x 2.0 x 0.75 inches
Casing	Black ABS Plastic / Aluminum
Warranty	1 Year