

HEMOSURE[®] iFOBTEST

One-step Immunological Fecal Occult Blood Test



“...in comparison with guaiac-based tests...
immunochemical tests are more patient-friendly...”
-American Cancer Society Guidelines for the Early Detection of Cancer, 2003

**Life just got easier
because
you deserve the best.**

Summary & Comparison

	HEMOSURE [®] FDA CLEARED & CLIA WAIVED	Traditional GUAIAIC TEST
Test Method	Immunoassay Test	Chemical Test
Interpretation	Clear accurate reading comparable to Pregnancy test	Blue color, often difficult to interpret
Lower GI Specificity	> 96%. Highly specific for human hemoglobin	Poor. False positivity from food and medicines
Overall Sensitivity	> 87%. Detects as low as 0.05ugHB/ml. specifically of human hemoglobin	50%. Detects 90 ugHB/ml. or higher of nonspecific hemoglobin
Accuracy	>94%	<86%
Dietary or Drug Restrictions or Preparations	NONE	YES. Several days.
Samples Required	ONE required, TWO recommended for CRC Screening*	THREE required
Reimbursement	Diagnostic code 82274QW: \$23.22 Screening code G0328QW:	Diagnostic code 82270: \$4.54 Screening code G0394:

*American Cancer Society, Can Colorectal Polyps and Cancer Be Found Early? 2/2007



HEMOSURE®

One-step Immunological Fecal Occult Blood Test

Product Numbers	T1-CK30 ***	T1-CK30T **	T1-CK50 **	PREPACK-iFOB20 ***	DUO-iFOB 20 *
Test Cassettes	T1-CT30 30	T1-CT30 30	T1-CT50 50	T1-CT20 20	T1-CT20 20
Collection Tubes	T1-TB50 50	T1-TB30 30	T1-TB50 50	● (Packed in Mailer)	● (Packed in Mailer)
Mailers	T1-CM50 50	●	●	20 (1 tubes per mailer)	DUO-CM10 10 (2 tubes per mailer)
	Best Value! 20 Tubes & 20 Mailers at No extra cost!		*** Complete Kit Ideal for Doctor's office ** Ideal for Laboratory use * DUO Package recommended for 2 consecutive samples.		

Hemosure Complete Product Offering:

First Sign™ Drug of Abuse Test	WH ACCU Test™ hCG Pregnancy Test	WH Ovultell™ LH Ovulation Test
Dipcard, cassette and built in cup format for detection of multiple drugs in urine	CLIA waived, cassette and dipcard format for detection of hCG in urine	Cassette and dipcard format for detection of LH in urine

To order product, contact your medical supply distributor. For more information on Hemosure® iFOB Test, call customer service at:

Ph: 318-227-4994 Fax 318-227-4996
www.americanscreeningcorp.com



HEMOSURE®

To Order Please Contact:
American Screening Corporation
 6658 Youree Dr., Ste 180 PMB 404
 "Complete Drug Testing Solutions"
 Shreveport, LA 71105
 Ph: 318-227-4994 Fax 318-227-4996
www.americanscreeningcorp.com



American Screening Corporation