



Clearly different.  
**Clearly better.**



Women's Health | Infectious Disease | Cardiovascular | Drugs of Abuse | GI Disorders

## Clearview® Strep A Exact II

### **Clearview Strep A Exact II**

Rapid technology provides direct detection of Group A Streptococcal Antigen.



### **Clearly different**

*Clearview's patented innovative technology provides highly sensitive and specific tests for Group A Strep.*

- **Rapid Response:** Positive result in as soon as 1 minute.
- **Simplicity:** CLIA waived, easy-to-use dipstick test.
- **Reliability:** Built-in control ensures accuracy.
- **Flexibility:** Results while patient waits.

### **Clearly better**

*Clearview Strep A Exact II provides physicians with fast, accurate results for an early diagnosis than can result in:*

- Shorter course of illness.
- Reduced severity of symptoms.
- Reduced spread of bacteria to others.
- Reduced likelihood of complications.




**Clearview Strep A Exact II**  
30 Dipsticks  
Product No. 92002  
CLIA Waived

**Clearview Strep A Exact II**  
30 Cassettes  
Product No. 92003  
CLIA Non-waived


## Clearview Strep A Comparison Grid

Product Name	Clearview Strep A Exact II Dipstick	Quidel QuickVue® Dipstick Strep A	Beckman Coulter ICON® Strep A	Genzyme OSOM® Ultra Strep A
Format	30 Dipsticks	50 Dipsticks	25 Dipsticks	25 Dipsticks
CLIA Status	Waived	Waived	Waived	Waived
Sample	Swab	Swab	Swab	Swab
Sensitivity	97%	92%	97%	92.6%
Specificity	95%	98%	95%	92.8%
Procedural Steps	5	5	5	10
Time to Result	5 minutes	5 minutes	5 minutes	7 minutes
Shelf Life	24 months	18 months	12 months	12 months
Storage	Room Temperature	Room Temperature	Room Temperature	Room Temperature

### Clearview Strep A Exact II Dipstick Procedure



### Clearview Strep A Exact II Cassette Procedure



For more information, visit:  
[www.americanscreeningcorp.com](http://www.americanscreeningcorp.com)



To Order Please Call:  
American Screening Corporation  
6658 Youree Dr Ste 180 PMB 404  
Shreveport, LA 71105  
Ph: 318-227-4994 \* Fax:318-227-4996  
Toll Free: 1-866-526-2873