

QUICKVUE[®]

One-Step hCG
Urine and Combo TEST

Two Simple Tests You Can Depend On For The Early Detection of Pregnancy



25 Test Urine Kit



50 Test Combo Kit

Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Strong positive may become visible sooner.

Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 25 mIU/mL sensitivity in serum or urine.

Convenient

- Room temperature storage.
- 18-month shelf life.

Easy to Use

- Simple, one-step procedure uses only 3 drops of sample.
- Two-color endpoint for easy-to-read and interpret results.


Reliable

- Highly sensitive assay with built-in controls.


Lot-to-Lot Consistency

- Highly controlled and monitored manufacturing environment.


americanpregnancy.com



Combo Procedure:



Urine Procedure:



To Order Please Contact:
American Screening Corporation
 6658 Youree Dr., Ste 180 PMB 404
 "Complete Drug Testing Solutions"
 Shreveport, LA 71105
 Ph: 318-227-4994 Fax 318-227-4996
www.americanscreeningcorp.com

Distributor Information: