

QUICKVUE[®]

Dipstick Strep A TEST

Test and Treat Today!

Group A Streptococcus, or Group A Strep, is a bacterium commonly found in the throat or on the skin. The bacteria can cause a range of infections from relatively mild sore throats and skin infections, to life threatening invasive disease. Group A infections can usually be treated with antibiotics and are spread by direct person to person contact.

The QuickVue Dipstick Strep A test is intended for the rapid qualitative detection of Group A Streptococcal

antigen from throat swabs or confirmation of presumptive Group A Streptococcal colonies recovered from culture.

The Dipstick format allows for an "off-line" extraction while maintaining its simplistic procedure, allowing the test to be granted CLIA-waiver. The two-colored endpoint and accuracy at the low levels of colony count make the QuickVue Dipstick Strep A test one of the most accurate and cost effective tests on the market today.



Rapid Diagnosis and Treatment

- Results in 5 minutes, or less.
- Patients can be tested and treated in the same office visit.
- CLIA waived

Easy to Use

- Dipstick format.
- Two-color endpoint for easy to read and interpret results.

Excellent Performance


- 92% sensitivity; 98% specificity
- 95% overall agreement.
- High sensitivity at low levels of colony count..

Convenient

- All components necessary to perform test are contained in kit.
- Room temperature storage (18 months from date of manufacture).
- Built-in Positive and Negative controls confirm the test has performed properly.
- Positive and Negative external controls provided in the kit.
- Individually packaged for convenient single-use testing and disposal.



Procedure:



Refer to package insert for complete instructions.



QuickVue Dipstick Strep A – 50 Test Kit:
Catalog #20108



QuickVue Dipstick Strep A – 25 Test Kit
Catalog #20125 (for sale outside the
U.S. only)



To Order Please Contact:
American Screening Corporation
6658 Youree Dr., Ste 180 PMB 404
“Complete Drug Testing Solutions”
Shreveport, LA 71105
Ph: 318-227-4994 Fax 318-227-4996
www.americanscreeningcorp.com