

## The #1 Selling PT/INR Monitor for Three Years in a Row.<sup>1</sup>

The INRatio2 PT/INR testing system combines reliable results with practical convenience making it an excellent in-office solution for anticoagulation patient management.



### On Board Quality Control...

Rigorous quality control checks are performed automatically with every test strip, ensuring a reliable result.



### Individually Wrapped Strips/ No Refrigeration Required...

Individually foil-wrapped strips can be stored at room temperature for up to 12 months.



### Small Sample Size...

1 unmeasured drop of blood applied directly to the strip performs both the test and quality controls.

## Reimbursement Information\*

| CPT Code    | Description  |
|-------------|--|
| 85610QW     | Prothrombin Time<br>(QW modifier indicates a CLIA waived test) |
| 99201-99205 | E&M New Patients   |
| 99211-99215 | E&M Established Patients                                       |
| 99363       | Anticoagulation Management<br>(initial 90 days of therapy)     |
| 99364       | Anticoagulation Management<br>(subsequent 90 days of therapy)  |

**\*NOTE:** Reimbursement may vary depending on your location/plan. It is the responsibility of the health care provider to properly code, document treatment performed, and seek reimbursement.

### References

- <sup>1</sup> HPIS comparative sales data, point of care PT/INR monitors sold Q4 2004 - Q2 2008.
- <sup>2</sup> The Seventh ACCP Conference on Anticoagulation Therapy: Evidence-Based Guidelines, Supplement to Chest, Volume 126, No. 3 (Supplement), September 2004: 2085-2095.
- <sup>3</sup> Chiquette E., Amato M., Bussey H., Comparison of an Anticoagulation Clinic with Usual Medical Care. Arch Intern Med. 1998; 158: 1611-1647

# For Ordering



6658 Youree Dr Ste 180 PMB 404  
Shreveport, La 71105

www.americanscreeningcorp.com  
Sales@americanscreeningcorp.com

Toll Free: 1-866-526-2873

Ph: 318-798-3386 Fax: 318-798-3306



Actual size



## True QC You Can See. Results You Can Trust.

Only the INRatio2 performs normal (QC1) and therapeutic (QC2) assays automatically on every test and provides quantitative, recordable results on the display—true QC you can see, start to finish.

### More Than Just A Humidity Check

These two levels of control use a true clotting control measure, following the standards of a split sample laboratory QC system, to provide assurance of testing integrity on every test.

### Easy Test Procedure

A simple fingerstick and 1 unmeasured drop of blood delivers results in about 60 seconds.



### Proven Adherence to Clinical Practice Guidelines

INRatio has always surpassed the ACCP guidelines for low ISI thromboplastin.<sup>2</sup>

### NEW and IMPROVED in INRatio2:

- Larger LCD screen and easy to read digits
- Simplified, icon based interface
- All kits include a complimentary tote
- Extended battery life supports two hundred tests
- Expanded test strip area for easier sample application
- Compatible with INRatio Test Strips
- Auto turn-on with strip insertion
- 33% smaller than INRatio

CLIA Waived

*INRatio2 will enable you to minimize the challenges of managing your oral anticoagulation patients.*



## INRatio2 Supports Your Goals of Improving Patient Care and Practice Efficiencies.

### Improve patient satisfaction...

Send patients home instead of to the lab for a PT/INR test.

### Reduce unnecessary delays...

Results in about 60 seconds instead of hours or even days.

### Allows prompt medication changes...

Counsel your patients face-to-face with dosage changes for improved TTR (Time in Therapeutic Range).<sup>3</sup>

### Maximize practice efficiencies...

Reduce time-consuming "call-backs", freeing your staff for more productive duties.

### E-Trainer...

A self-paced, interactive, web or CD based educational tool. Learn more at [www.hemosense.com](http://www.hemosense.com) today!

#### INRatio2 Pro Kit Cat#0200431

- |                         |                              |
|-------------------------|------------------------------|
| 1 INRatio2 Monitor      | 1 Training Video             |
| 1 User's Guide          | 10 Blood Lancets             |
| 1 Quick Reference Guide | 10 MicroSafe Capillary Tubes |
| 1 Carrying Tote         |                              |
| 1 Power Supply          |                              |

#### Test Strip Bulk Pack

Cat#0100139  
48 INRatio Test Strips Per Box