


# OraLine®

The **ONLY** instant oral fluid drug screen that meets the SAMHSA cutoff for THC\*

## THC 4ng/mL (Marijuana)



No additional oral fluid needed for confirmation

### Oraline Cutoff

THC 4 ng/mL  
COC 25 ng/mL  
MET 50ng/mL  
OPI 40ng/mL

1



4 green lines indicate positive readings.

2



After specimen collection, a red line opposite the control "C" symbol indicates the screen is valid.

3



Red lines of any color indicate a negative result (no drug detected).

4



Example above: The screen is negative for opiates and methamphetamine but positive for THC and cocaine. Total test time is 10 minutes.

