

































← **FEED DIRECTION**

	TEST AND READING TIME	ABNORMAL (LOW)	NORMAL	ABNORMAL (HIGH)
	Creatinine 45 seconds	 Negative  10	 20  50  100  200 mg/dl	
	Nitrite 45 seconds		 0  0.1-0.2  0.5-5.0	 >15mg/dl
	Glutaraldehyde 45 seconds		 Negative	 Positive
	pH Immediate	 2  3	 4  5  7  9	 ≥10
	Specific Gravity 45 seconds	 1.000	 1.005  1.015  1.025	 ≥1.030
	Bleach 30 seconds		 Negative	  Positive 
	Pyridinium Chlorochromate 30 seconds		 Negative	  Positive 
			H A N D L E	

LTX Item Number 000629-001
ZPL1-DA6